

# Cav902B

## LIPOLASER

### User&Training Manual



# 902B

## LIPOLASER

### User & Training Manual

---

## Diode Laser Safety Guidelines

---

The Following information is provided for the correct utilization of IPL Skin treatment system. The information includes not only the accident protection regulations the products comply with, but also the effective regarding proper use of the products.

### 1) Electric safety regulation

The EMC Directive 89/336/EEC amended by 92/31/EC & 93/68/EEC EN 55015:2000+A1+A2/LUM200704025/EMC EN 61000-3-2:2000+A2/LUM200704025/EMC EN 61000-3-3:1995+A1+A2/LUM200704025/EMC EN 61547:1997+A1/LUM200704025/EMC

### 2) Material ROHS test 2002/95/EC

Restriction of use of certain Hazardous

Substance in Electrical & electronics Equipment/ROHS.

3) Electromagnetic radiation safety regulation LVD/EN 60601-1:2006/AGC10430904GZ05E3EMC/EN 60601-1-2:2001+A1:2006/AGC10430904GZ05E1

4) ISO 13485:2003 MEDICAL/33021200005459

5) ITC MEDICAL CE

IPL SKIN TREATMENT SYSTEM belongs to II a according to MDD directive, It designed according to EN 980:2003; EN 1041:1998; EN ISO 14971:2007; EN 60601-1-2:2007; EN 61000-3-2:2006; EN 61000-3-3:1995+A1:2001+A2:2006; EN 60601-2-22:1996; EN 60601-1-4-1996+A1-1999; EN 60601-1:2006.

**902B**

**LIPOLASER**

**User & Training Manual**

---

## **1. LIPOLASER common sense**

---

**Is LipoLaser Safe?**

Yes, LipoLaser treatments are completely non-invasive and clients feel nothing. LipoLaser's low level laser technology has been studied and used for over 30 years for many biomedical purposes. FDA approved for pain and inflammation, doctors can use it for fat loss .

**Do I have to do other things to lose the fat?**

Only for bigger results- 1 inches without changes, 5 inches with major lifestyle changes, like diet and exercise - we have support if you choose to change your life while doing the sessions

**What do you do to keep the fat off?**

Normal things you need to keep weight off - diet and exercise, drinking water, etc.

**Where Does the Fat Go?**

That fat can be burned as energy in the body, or normal metabolic processing—just like when you exercise and release fat

**How Quickly Will I See Results**

Many clients see measureable results after their first session.

**What Kind of Results are Produced?**

1-5 inches average-the most lost in 9 sessions is 23 inches  
People who make minor lifestyle changes, are on the higher end of the average, and those who don't, the lower end.

**How many times do I need it ?**

3 times per week, 9 times per course of treatment.

# 902B

## LIPOLASER

### User&Training Manual

---

## 2.Product Information

---

### 1.PRODUCT INFORMATION

#### APPLICATION

- \*Intensive physical lipolysis to remove fat
- \*Excess fat cell melted
- \*Body slimming, cellulite reduction
- \*Smooth fatigue
- \*Remove obstruction from channels and collaterals
- \*Promote and accelerate the body's metabolism



#### FEATURES

- \*Laser wave length:635-650nm
- \*Maximum output power:400W
- \*Single Output:100mW
- \*Numbers of Treatment Probe:8
- \*Numbers of Diode laser:62
- \*Timer:30minutes MAX
- \*Input current:240v or 110V AC
- \*Working Temperature:10-40degree
- \*Weight:21Kg
- \*Package size:58\*50\*33cm
- \*Certifications:CE LVD,CE EMC


# 902B

## LIPOLASER

### User&Training Manual

### 3.Installation Alignment and operation procedures


3.1 Open the Package and connect all parts.



This is 2pinsocket,connect to ellipse Pad

This is 4pinsocket,connect to rectangle Pad

6Piece Laser Pad with 10lamp



2Piece Laser Pad with 2lamp




# 902B

## LIPOLASER

### User&Training Manual

### 3.Installation Alignment and operation procedures



3.1. Turn the key on the key switch and urgency Button to power on the machine, and the led on the screen will display.

3.1.Adjust the timer for big probe,then press On to make 6 piece big probe working.

3.2.Press the on/off button to control the small probe working.

# 902B LIPOLASER


## User&Training Manual

---

### 4. Treatment procedure

---

4.1 Cleanse the skin, Put the Lipolaser probe to body where you want lost fat. then use a belt keep the pad to body.





4.2 You relax in a comfortable position while a trained LipoLaser technician places two multi-laser paddles and two smaller lasers to the target areas. Then lay back, relax, read, or whatever you like.

---

### 5. Before and after compare

---



# 902B

## LIPOLASER

### User&Training Manual

---

## 6.0ur Diode Laser Device superiority

---

6.1.Whole ABS Plastic Mould.

6.2.It does not present any contagious risk

6 .3.stable laser energy output

6.4.Good energy and effect

6.5.Super treatment effect

**We Are Looking for Sale agent over the world.**

